QUICK START GUIDE ECLIPSE 5°



Unpack Your Eclipse 5



Wall (AC) Power Supply

Eclipse Cart



Motor Vehicle (DC) Power Supply




Battery (Power Cartridge)



CAIRE[®] Call Us Toll-Free Today! **1.833.203.2473** www.cairemedical.com

QUICK START GUIDE ECLIPSE 5"



1.833.203.2473 www.cairemedical.com



Do not smoke near unit or while operating unit.

QUICK START GUIDE ECLIPSE 5"

-		0	1	D	h	4
0	н.	6	L_	Γ.	Ľ	1

Flow Setting Indicator: This is the main focus on your control panel. Your home care provider will correctly set your prescribed flow for the Continuous Flow Mode (LPM) and/or your Pulse Dose Mode (mL) settings. Each time you power the device ON, the previous mode and/or setting has been saved and will be used at start-up.

Power Cartridge (battery) Status Gauge: This indicator displays the charge remaining in the battery. Each of the five horizontal gray bars represents approximately 20% of the total battery charge. When the battery is being charged, the charge indicator bars will blink in a waterfall-type fashion. If the battery is not installed, or if it is improperly installed, the Power Cartridge (battery) Status Gauge will not be illuminated or will flash on and off.

Buzzer: An audible alarm (or buzzer) alerts you to the operating condition of the device, either a warning or failure, and confirms a valid key press.



AIR

Powering On the Eclipse 5

To connect to AC power source: Insert the AC Power Supply cord securely into the External Power Receptacle on the right side of the device. If your power supply has three pins, plug the AC power cord into a grounded AC outlet. If your power supply has a two pin plug, plug the AC cord into a grounded or non-grounded AC outlet.





When the Eclipse 5 is plugged in properly, a green indicator on the AC Power Supply will light and the External Power Indicator, will appear on the user control panel.

Install the Battery: Slide the battery into the empty compartment until it locks into position and is flush with the back of the device.




When the battery is properly installed and the Eclipse 5 is ON, the Power Cartridge Status Gauge **will appear**.





ON/OFF Button (Green) Indicator


Press and hold the "ON/OFF" button for two (2) seconds to power ON your Eclipse 5. A brief audible and visual power on self-test will occur. Verify that all the indicators are illuminated and the buzzer sounds for three (3) seconds. Press and hold the "ON/OFF" button for two (2) seconds to power OFF your Eclipse 5.



Press the Flow Mode button to select your desired flow mode. Pressing this button repeatedly will toggle you back and forth between continuous flow and pulse dose modes. In Continuous Flow Mode the oxygen is provided at a constant flow rate between 0.5 and 3.0 LPM. In Pulse Dose Mode, oxygen is supplied in a bolus at the beginning of each inspiration, providing a selectable range of settings 1-9 with bolus sizes 16-192 mL. Flow settings in both modes should be prescribed for you by your physician.

Call Us Toll-Free Today! **1.833.203.2473** WWW.CAIREMEDICAL.COM

QUICK START GUIDE ECLIPSE 5°




Using the + or - buttons, adjust the setting to the prescription supplied by your physician.

Continuous Flow Mode Display



Pulse Dose Mode Display



WARNING: IT IS VERY IMPORTANT TO SELECT ONLY THE PRESCRIBED LEVEL OF OXYGEN. DO NOT CHANGE THE FLOW SELECTION UNLESS YOU HAVE BEEN DIRECTED TO DO SO BY A LICENSED CLINICIAN. THE PORTABLE OXYGEN CONCENTRA-TOR MAY BE USED DURING SLEEP UNDER THE RECOMMENDATION OF A OUALIFIED CLINICIAN.



-

Begin Using the Eclipse 5

Connect the oxygen supply tubing to the oxygen outlet and connect the oxygen cannula per the instructions provided with the cannula, or per the instructions from your home care provider.



Air Inlet Filter





Your unit is equipped with an Air Inlet Filter. Check daily and clean once a week. Breathe normally through the nasal cannula.

Universal Cart Operation

Align the locater holes on the bottom of the device to the tabs on the Universal Cart. Align the large threaded screw to the insert at the back of the Eclipse 5. Then hand-tighten the knob on the cart until the Universal Cart is tightly secured to the device. Depress the push-button on the telescoping Universal Cart handle to adjust the height of the handle.





THIS QUICK START GUIDE IS NOT INTENDED TO REPLACE THE DETAILED INFORMATION AVAILABLE IN YOUR USER MANUAL, SO PLEASE REFER TO THE USER MANUAL FOR IMPORTANT CRITICAL INFORMATION REGARDING USE AND OPERATION OF THE PRODUCT.

PLEASE CONSULT THE APPLICABLE PRODUCT INSTRUCTIONS FOR USE FOR PRODUCT INDICATIONS, CONTRAINDICATIONS, WARNINGS, PRECAUTIONS, AND DETAILED SAFETY INFORMATION.

© Copyright 2017 Chart Industries ML-CONC0062 F



WWW.CAIREMEDICAL.COM

AIR