



Presented With Compliments From
www.DirectHomeMedical.com

AirFit™ P10... a quieter mask to fit your sleep.

The ResMed AirFit™ P10 nasal pillows mask delivers the premium performance and reliability of the Swift™ FX mask, yet it's approximately 50% lighter and 50% quieter.¹ A recent study found that patients slept more than 40 minutes longer per night using the AirFit P10 compared to the leading nasal pillows mask.²



Benefits for you ...

Supports patient acceptance thanks to ease of use, lighter weight and greater comfort.

QuickFit™ auto-adjusting headgear can be fitted in 30 seconds or less.³

Simple sizing and color-coded pillows make fitting easier, encouraging patient compliance.

... and your patients

QuietAir™ woven-mesh vent makes this ResMed's quietest mask at just 21.4 dBA.*²

Minimal design and light weight (1.6 oz) deliver a high level of comfort and facial freedom.

Just three parts – frame, cushion and headgear – make it easy to assemble and take apart for cleaning.


* Declared dual numbered noise emission values (with uncertainty 3dBA) in accordance with ISO 4871.

¹ ResMed Ref: D000-0201; eLB0798-206.

² ResMed Clinical Study FECS3 Ref: D629-117 (page 11). ResMed objective data (pressure, leak, usage, AHI).

³ ResMed Clinical Study FECS3 Ref: D629-117 (page 9). 90% of participants fitted the mask comfortably in less than 30 seconds.

AirFit™ P10
Nasal Pillows System



ResMed

Comfort with stability.

A recent clinical study found that more than 70% of patients preferred the new ResMed AirFit P10 to the popular Swift FX mask for both comfort and stability. On average, participants also slept more than 40 minutes longer per night using AirFit P10 compared to the leading nasal pillows mask.²

The AirFit P10 nasal pillows mask provides a stable seal and a high level of comfort, with excellent usability. It also provides patients with a real sense of facial freedom. The mask's innovative QuickFit elastic headgear is unobtrusive yet provides a two-point head support with minimal facial contact—designed to encourage easy patient acceptance.



Flexible pillows

with trampoline action.

Dual-wall nasal pillows

seal on contact and offer a wide range of sizes.

Color-coded cushions

for easy identification.

Three-part simplicity

for fast and easy assembly.

Lightweight integrated flexi-tube allows for flexibility without sacrificing stability.

QuickFit elastic headgear

is auto-adjusting and can be fitted in 30 seconds or less.²

QuietAir woven-mesh venting gently directs air away from you and your partner.

~50% lighter than the popular Swift FX mask