

Enjoy comfort *and* a great seal with the world's first nasal mask featuring both Custom Fit Technology and an Active Headgear Connector. With three sizes and a Flexible Spacebar, Aura will give you the rest you deserve.

## Enter a New Aura of Comfort *and* Seal

### AIR<sup>°</sup>gel

Maximum comfort, minimum weight. Our patented AIR<sup>°</sup>gel minimizes pressure points and skin irritation.



### Quick-Release Headgear

Lightweight, breathable headgear slides off easily with one hand.

### Custom Fit Technology<sup>™</sup>

The shapeable shell of each mask can be easily molded to fit individual needs. Flex, fit, forget.

### 360° Rotating Swivel

For optimal tube positioning and easier mobility.

### Flexible Spacebar<sup>™</sup>

Unique bendable design provides superior stability, holding a seal without contacting the skin.

### Size Indicator

available in three sizes. **S/M/L**

### Advanced Cushion Technology<sup>™</sup> (A.C.T.)

Our ergonomic cushion technology delivers an effective seal for the most comfortable fit.

### Active Headgear Connector<sup>™</sup>

Moves as you move for a reliable seal.



Made in the USA  
of US and imported parts

## What makes our Aura Nasal Mask different?

- Custom Fit Technology is unique to Sleepnet. Our soft-shell mask has a shapeable wire allowing for optimum flexibility, fit, and comfort.
- Our Active Headgear Connector moves when you move to help maintain a seal throughout the night.
- Our AIR<sup>o</sup>gel with Advanced Cushion Technology creates an ergonomic gel cushion that reduces pressure points on the skin.
- Available in three sizes with simple and breathable quick-release EZ Fit headgear.

### About Advanced Cushion Technology (A.C.T.)

Our patented AIR<sup>o</sup>gel with an ergonomic streamlined design and luxurious feel reduces pressure while maintaining an effective seal.



### How does Aura compare to the competition?

Features	Sleepnet Aura Mask	ResMed Mirage™ FX Mask	Fisher & Paykel Eson™ Mask	Philips Respironics Wisp™ Mask
Custom Fit Technology	●			
Advanced Cushion Technology	●			
No Assembly Required	●			
Active Headgear Connector	●			
Patented AIR <sup>o</sup> gel	●			
Breathable Headgear	●	●	●	●

### Aura Mask with EZ Fit Headgear Specifications:

Product Number	Style	Size	Mask Weight	Mask Dimensions
50176	Vented	S	62 g	5.50"D x 2.85"W x 3.80"H
50177	Vented	M	68 g	5.50"D x 3.20"W x 4.00"H
50178	Vented	L	72 g	5.50"D x 3.45"W x 4.00"H
55051	EZ Fit	One Size	N/A	N/A



This product is not manufactured with natural rubber latex.

CE marking expected late 2013  
Pending FDA Clearance

Also from Sleepnet:

## Mojo<sup>®</sup>

Your Mask, Your Way

Our full face mask also features AIR<sup>o</sup>gel.



Contact us today.

To place your order, call an authorized distributor.



Presented With Compliments From  
[www.DirectHomeMedical.com](http://www.DirectHomeMedical.com)