

Nasal-Aire® II CPAP Interface

Something New is in the Aire...

The Nasal-Aire II is the **next generation of the ever-popular Nasal-Aire** interface. The Nasal-Aire II shares many of the **comfort and compliance benefits** of the original Nasal-Aire, adding improved performance and a reduction in sound level.

This lightweight interface is comfortable for patients, simple to use and easy to care for. Nasal-Aire II is **sold with all six sizes for easy fitting** and includes headgear for stability.



Nasal-Aire II
CPAP Interface
(K2A)

Performance Features

- Ergonomically designed to fit any face with virtually no leaks
- Virtually eliminates uncomfortable pressure points on the face allowing the interface to be worn for any duration without discomfort
- Secure fit design **allows the user to sleep in any position**
- Minimal design **reduces feelings of claustrophobia**
- New ultra-soft nasal prong material for improved patient comfort

Design Features

- 100% latex free
- Medium or large size fits ninety percent of patients
- Can be worn over the ears or down the chest
- Small and quiet with virtually no noise
- No nasal bridge or forehead contact **eliminates nasal bridge sores and air leaks near eyes**
- Swivel tubing coupler to accommodate free movement and optimal performance during the night
- Supplemental oxygen ports

Reimbursement Information

Product Description	Code Description	HCPCS Coding	Frequency	Allowable Range
NASAL-AIRE II INSERT	Cushion for nasal application device	A7032	Billing at initial setup Two per month	\$34 -\$40
NASAL-AIRE II APPLICATION	Nasal interface (mask or cannula type) used with positive airway pressure device, with or without head strap	A7034	Billing at initial setup One per 3 months	\$99 -\$117
HEADGEAR	Headgear used with positive airway pressure device	A7035	Billing at initial setup One per 6 months	\$33 -\$39

Accessories and Replacement Parts

PART NUMBER	PRODUCT DESCRIPTION	PART NUMBER	PRODUCT DESCRIPTION
K2A	Nasal-Aire II All Size Kit, Includes XL, LG, MP, MD, SM, XS*	LG05	Nasal-Aire II Cannula Only, LG
K2LG	Nasal-Aire II New Patient Setup Kit, MP, LG, XL*	MD03	Nasal-Aire II Cannula Only, MD
K2MD	Nasal-Aire II New Patient Setup Kit, SM, MD, MP*	MP04	Nasal-Aire II Cannula Only, MP
K2MP	Nasal-Aire II New Patient Setup Kit, MD, MP, LG*	SM02	Nasal-Aire II Cannula Only, SM
K2SM	Nasal-Aire II New Patient Setup Kit, XSM, SM, MD*	XL06	Nasal-Aire II Cannula Only, XL
PA401	Nasal-Aire II Petite, Size A	XS01	Nasal-Aire II Cannula Only, XS
PB402	Nasal-Aire II Petite, Size B	NAT320	Nasal-Aire II Frame Kit
PC403	Nasal-Aire II Petite, Size C	X2105	Nasal-Aire II Headgear with OHS
PD404	Nasal-Aire II Petite, Size D	Assembly	
PE405A	Nasal-Aire II Petite, Size E	P430	Nasal-Aire II Petite – Headgear
LG305	Nasal-Aire II, LG*	PA411	Nasal-Aire II Petite – Cannula Only, Size A
MD303	Nasal-Aire II, MD*	PB412	Nasal-Aire II Petite – Cannula Only, Size B
MP304	Nasal-Aire II, MP*	PC413	Nasal-Aire II Petite – Cannula Only, Size C
SM302	Nasal-Aire II, SM*	PD414	Nasal-Aire II Petite – Cannula Only, Size D
XL306A	Nasal-Aire II, XL*	PE415A	Nasal-Aire II Petite – Cannula Only, Size E
XS301	Nasal-Aire II, XS*	X2007	Sizing Gauge – Nasal-Aire II
		X4007	Sizing Gauge – Nasal-Aire II Petite

* Includes headgear, frame, bolo and multiple nasal cannulas.

** Includes headgear, frame, nasal cannula and bolo.



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com