

Mask and Accessories Materials List

All Parts are LATEX FREE

Product	Material
ComfortGel™	
Outer Flap	Silicone/Polycarbonate
Gel Cushion Outer Shell	Polyurethane Film
Cushion Gel	Polyurethane Gel/EVA
Swivel Elbow	Polycarbonate
Faceplate	Polycarbonate
Elbow Split Washer	Polycarbonate and/or PETG
Forehead Support Arm and Bracket	Polycarbonate
Forehead Cushion	Silicone
Swivel Headgear Clip	Acetal
Premium Headgear	UBL/HyPUR-cel® foam/Lycra®

ComfortFull™ 2	
Cushion	Silicone
Faceplate	Polycarbonate
Valve Body	Polycarbonate
Forehead Cushion	Silicone
Premium Headgear	UBL/HyPUR-cel® foam/Lycra®

ComfortLite™ 2	
Cushion (Simple or DirectSeal)	Silicone
Simple Cushion Elbow	Polycarbonate
DirectSeal Cushion Cradle	Polycarbonate
DirectSeal Cushion Support Band	Aluminum
Crown Swivel	Polycarbonate
Flex Tube	EVA copolymer
Pivot Tube Insert	Polycarbonate
Height Adjustment Nut	Acetal
Height Adjustment Tube	Polycarbonate
Pivot Mount Strap	Polycarbonate
Pivot Tube	Polypropylene
Headgear Halo	Nylon
Headgear Padding	Nylon/Lycra® Laminate
Headgear Strap	Nylon/Lycra® Laminate

ComfortSelect™	
Cushion	Silicone
Split Washer	Polycarbonate and/or PETG
Swivel Elbow	Polycarbonate
Faceplate	Polycarbonate
Forehead Cushion	Silicone
Deluxe Headgear	Velstretch®/Lycra® Laminate

Product	Material
ComfortLite™	
Cushion (Simple or DirectSeal)	Silicone
Simple Cushion Elbow	Polycarbonate
DirectSeal Cushion Cradle	Polycarbonate
DirectSeal Cushion Support Band	Aluminum
Crown Swivel	Polycarbonate
Flex Tube	EVA copolymer
Pivot Tube Insert	Polycarbonate
Height Adjustment Nut	Acetal
Height Adjustment Tube	Polycarbonate
Pivot Mount Strap	Polycarbonate
Pivot Tube	Polypropylene
Headband and Crown Strap	High Density Polyethylene
Headband and Crown Strap Liner	Coolmax Mesh, Polyester Foam, Tricot
Rear Strap	Velcro® Hook and Loop
Rear Strap Clip	High Density Polyethylene

ComfortCurve™	
Cheek Pads	Silicone
Cheek Pads	Polyurethane Foam
Cushions	Silicone
Bridge	Thermoplastic Polyester Elastomer
Cheek Pad Supports	Polycarbonate
Headgear Strap	Velcro® Velstretch® Nylon Loop
Headgear Clip	Polycarbonate
Elbow	Polycarbonate

ComfortClassic™	
Cushion	Silicone
Faceplate	Polycarbonate
Split Washer	Polycarbonate and/or PETG
White Port Cap	Flexible Vinyl (PVC DEHP Free)
Swivel Elbow w/Exhalation Port	Polycarbonate
GEL Forehead Spacer	Polyurethane Elastic Film & GEL
Deluxe Headgear	Velstretch®/Lycra® Laminate

Product	Material
Profile™ Lite/GoldSeal®	
Cushion Outer Shell	Polyurethane Elastic Film
Cushion Gel	Polyurethane Gel
White Port Cap	Flexible Vinyl (PVC)
Split Elbow	Polycarbonate and/or PETG
Swivel Elbow	Polycarbonate
Faceplate	Polycarbonate
Spacer	White Polyethyl Foam
Deluxe Headgear	Velstretch®/Lycra® Laminate

Simplicity™ Nasal Mask	
Faceplate	Polycarbonate
Tubing	Ethylene Vinyl Acetate (EVA)
Cushion	Silicone
Simplicity Headgear	Neoprene/Nylon
Deluxe Softcap	Velstretch®/Lycra® Nylon

Total™ Face Mask	
Cushion	Silicone
Faceplate	Polycarbonate
Seal Ring	Thermoplastic Polyester Elastomer
Valve Body	Polycarbonate
Total Headgear	Neoprene/Nylon

Spectrum® Reusable Full Face	
Cushion	Silicone
Cushion	Polyurethane
Faceplate	Propionate
Seal Ring	Kraton Elastomer
Valve Body	Polycarbonate
Forehead Pad	White Polyethyl Foam
Comfort Flap	Polyurethane

Product	Material
Headgear	
Deluxe Headgear	Velstretch®/Lycra® Laminate
ComfortHeadgear™	Nylon/Lycra® Laminate
Premium Headgear	UBL/HyPUR-cel® foam/Lycra®
Simple Strap	
Mesh	Polyester
Strap Loop/Hook	Velstretch®/Nylon Hook
Softcaps	
Cap	Lycra®
Strap Loop/Hook	Velstretch®/Nylon Hook


RESPIRONICS®

Customer Service: 800-345-6443 or 724-387-4000
Respironics Europe, Africa, Middle East: +33-1-47-52-30-00
Respironics Asia Pacific: +852-3194-2280

©2006 Respironics, Inc. and its affiliates. All rights reserved.

Respironics is a registered trademark; Comfort Series, ComfortCurve, ComfortLite, ComfortGel, ComfortSelect, ComfortClassic, ComfortFull, ComfortHeadgear, Softcap, SimpleStrap, Profile Lite, Simplicity, Total, Spectrum, and Contour Deluxe are trademarks of Respironics, Inc. and its affiliates. Velstretch is a trademark of Velcro Industries. Lycra is a registered trademark of INVISTA Company. Akwadyne is a registered trademark of Comfort Technologies, Inc. HyPUR-cel is a registered trademark of Rubberlite, Inc.