



THE NEW EVERFLO™ Q IS ULTRA QUIET

NOT TO MENTION LIGHT, SMALL AND REMARKABLY AFFORDABLE

Respironics' compact, 31-pound EverFlo Q stationary oxygen concentrator is the first ultra-quiet concentrator available at an affordable price. The new EverFlo Q combines the small size, light weight and low maintenance of our popular EverFlo oxygen concentrator, with the added value of an ultra-quiet sound level. With EverFlo Q, you get the benefits of an ultra-quiet concentrator without the burden of paying an ultra-large premium. And, that's good for you and your patients.



Presented With Compliments From
www.DirectHomeMedical.com



EverFlo™ Q
5 LPM Stationary Concentrator

EverFlo™ Q

5 LPM Stationary Concentrator

With a sound level of 40 decibels, the new EverFlo Q can significantly reduce noise complaints while giving you a feature advantage that is attractive to patients and physicians. EverFlo Q is small and light – translating into lower shipping and storage costs and reduced risk of injury. Because EverFlo Q has only two filters, you can spend less time and money on regular maintenance or repair calls.

The easy-to-use EverFlo Q is also attractive to your patients because it produces less heat and requires less electricity to operate. So, now you can control oxygen therapy costs and deliver the features you and your patients want.



PRODUCT SPECIFICATIONS

Sound Level	40 dBA typical
Input Voltage	120 VAC (+/- 10%)
Input Frequency	60 Hz
Average Power Consumption	350 watts
Oxygen Concentration* (at 5 LPM)	93% (+/- 3%)
Liter Flow	0.5-5 liters per minute
Weight	31 lbs (14 kg)
Dimensions	23 in (584 mm) H, 15 in (381 mm) W, 9.5 in (241 mm) D
OPI (Oxygen Percentage Indicator) Alarm Levels	Low Oxygen: 82%, Very Low Oxygen: 70%
Operating Temperature	55°F to 90°F (12°C to 32°C)
Storage/Transport Humidity	-30°F to 160°F (-34°C to 71°C) up to 95% relative humidity
Operating Humidity	Up to 95% relative humidity
Outlet Pressure	5.5 PSI
Operating Altitude	0 to 7,500 ft (0 to 2,286 m)

*Device operation above or outside of the specified voltage, LPM, temperature, humidity and/or altitude values may decrease oxygen concentration levels.

ORDERING INFORMATION

CONCENTRATOR SYSTEMS	PART NUMBER
EverFlo Q Oxygen Concentrator	1020015
EverFlo Q Oxygen Concentrator with OPI (Oxygen Percentage Indicator)	1020014

CONCENTRATOR ACCESSORIES

Compressor Intake Filter	1038831
Micro Disk Filter	520
Humidifier Connector Tube	1039642
Stand-alone Pediatric Flow Meter	H537-1
Integrated Pediatric Flow Meter	H644
Flow Meter Locking Kit	528
Caster Brake Kit	H649
EverFlo User Manual in English, French, Spanish and Portuguese	1041063
EverFlo Service Manual	1039055

EXTENDED WARRANTY

Two-Year Extended Warranty with Compressor	1038571
--	---------

KEY FEATURES:

- Ultra quiet – 40 dBA typical
- Premium sound level with affordable pricing
- 31 pounds, small and compact
- Easy to store and transport
- No provider filter change for two years
- Only two filters to maintain
- OPI and non-OPI versions

HOW THE EVERFLO Q CAN BENEFIT YOUR OXYGEN BUSINESS:

- Attract new patients with quiet, small stationary concentrator
- Increase patient satisfaction by providing best combination of sound, size and ease-of-use in a stationary concentrator
- Decrease the time and cost of addressing noise complaints
- Differentiate your oxygen business to physicians

RESPIRONICS®

Envisioning tomorrow. Improving today.™

Customer Service:
800-345-6443 or 724-387-4000
Respironics Europe, Asia, Middle East:
+33-1-47-52-30-00
Respironics Asia Pacific:
+852-3194-2280

www.respironics.com

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician.

Respironics is a registered trademark, EverFlo is a trademark and *Envisioning tomorrow. Improving today.* is a service mark of Respironics, Inc. and its affiliates.

©2008 Respironics, Inc. and its affiliates. All rights reserved.

Geyer WMB 12/21/07 MCI 4101504 PN 1050326