

INOGEN ONE[®] G2 System



Presented With Compliments From
www.DirectHomeMedical.com



Inogen One[®] G2
Universal Power Supply*
Item # BA107/207



Inogen One[®] G2
Cart & Carry Bag*
Item # CA200



Inogen One[®] G2
12-Cell Lithium Ion Battery*
Item #BA200



Inogen One[®] G2
24-Cell Lithium Ion Battery**
Item #BA224

The Inogen One[®] G2 System (IS-200):

- 40% Smaller
- 25% Lighter
- 20% More Oxygen
- Authorized by FAA onboard aircraft

*Included with system

**Accessory sold separately

Oxygen. Anytime. Anywhere®

The new Inogen One® G2 delivers the independence of our original Inogen One portable concentrator, but it's smaller, lighter, and makes more oxygen with a longer battery life. It's a portable concentrator small enough to hold in your hand. That means there are no tanks to refill and you never have to worry about running out of oxygen.

The Inogen One® G2 is your one solution, your single solution, for oxygen at home, away, and for travel (approved for use on all US airlines). The Inogen One® G2 is Oxygen. Anytime. Anywhere®.

To learn more, call toll-free 1-877-446-6436.

Weight	7.25 lbs. (includes 12-cell battery)
Size	Length: 10.7" Width: 3.9" Height: 9.5" (including handle)
Oxygen Flow	Pulse dose delivery system Five flow settings from 1-5
Power Sensing	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power: for mobile use in car/airplane
Battery	Duration (12-cell battery): 3.5 to 4 hours* Recharge approximately 4 hours with AC or DC power Duration (24-cell battery): 7 to 8 hours* Recharge approximately 8 hours with AC or DC power
Noise Level	38 Decibels*
Warranty	3 year limited warranty**
Operation	Simple control functions and easy-to-read LCD display
Use	Designed for stationary, portable, and travel

*At setting 2

**See Inogen One® G2 Limited Warranty Statement for complete warranty information

Manufactured by:



326 Bollay Drive
Goleta, California 93117
1-877-446-6436
www.inogenone.com
info@inogen.net



MKT-0238A