

DreamFit™

NASAL MASK



Presented With Compliments From
www.DirectHomeMedical.com

New Easy-to-Fit CPAP Interface

The new *DreamFit*™ Nasal Mask from Puritan Bennett expands our selection of CPAP interfaces. Now you can offer a cost-effective nasal mask option when equipping patients with Puritan Bennett® sleep therapy devices (see reverse side). The *DreamFit* Nasal Mask has a number of features to enhance comfort and accommodate a range of sleeping positions. It is designed with a flexible “bellows” chamber between the mask and interface that swivels to accommodate patient movement and reduce leaks. It also features a single closure, so it’s easy for patients to put on and take off. The soft, pliable *DreamFit* mask creates an effective, yet comfortable seal.



Puritan Bennett has pioneered innovative interface technologies, such as the Nasal Pillows used in our *Breeze SleepGear* CPAP Interface System.

Unique bellows chamber swivels to help prevent leaks during patient movement.

Soft Breath-O-Prene® headgear has simple, one-step closure.

Simple, open design minimizes contact to the face.

Like our popular *Breeze SleepGear*™ CPAP Interface System, the airway hose routes over the head to enable sleep in virtually any position.

Comes in three sizes (standard, large, shallow) to fit a wide range of patients.

DreamFit

NASAL MASK

Specifications

Materials of Construction

- *DreamFit* Mask: Perimeter seal is medical-grade silicone
- Straps: Breath-O-Prene
- All *DreamFit* components are latex-free

Pressure Range

- 3 – 35 cm H₂O

Hose Attachment

- Fits 22 mm O.D. connection

Ordering Information

Description	Catalog Number
<i>DreamFit</i> Interface with Large Mask	DF301
<i>DreamFit</i> Interface with Standard Mask	DF302
<i>DreamFit</i> Interface with Shallow Mask	DF303
<i>DreamSeal</i> Adapter Kit	10001990
<i>DreamFit</i> Strap Kit	10001994
<i>DreamFit</i> Hose Kit	10002002
<i>DreamFit</i> Silicone Parts Kit	10002003

World's Smallest, Lightest Sleep Therapy Devices

Let sleep apnea patients experience the benefits of the world's smallest, lightest CPAP device. Puritan Bennett *GoodKnight*® 420 Series devices are only 3" tall, 5.5" wide and 8" deep, making it easy for patients to comply with therapy wherever they go. These reliable units are whisper-quiet and have low air resistance to enhance patient comfort. The *GoodKnight* product family includes:

- ***GoodKnight* 420G CPAP System**
Entry-level model
- ***GoodKnight* 420S CPAP System**
Includes compliance tracking features
- ***GoodKnight* 420E Auto-CPAP System**
Auto-CPAP function automatically adapts to airway pressure needs
- **(New!) *GoodKnight* 425 Bi-Level® Device**
Delivers separate inspiratory and expiratory pressure levels
- **(New!) *GoodKnight* H₂O Heated Humidifier**
Features a small footprint that fits neatly with any *GoodKnight* 420/425 sleep therapy device



tyco

Healthcare

**Puritan
Bennett**

4280 Hacienda Drive
Pleasanton, CA 94588
Tel 925.463.4000
Toll Free 1.800.635.5267
www.puritanbennett.com

Tyco Healthcare UK LTD.
154 Fareham Road
Gosport, UK PO13 0AS
Tel +44.1329.224000

Caution: Federal (U.S.A.)
law restricts this device
to sale by or on the order
of a physician.

Specifications subject to
change without notice.

© 2005 Nellcor Puritan
Bennett Inc. All rights
reserved.

Breath-O-Prene is a
registered trademark of
AccuMED Technologies, Inc.

A.ae 2175-0305
ST06604