


Interfaces

FlexiFit™ 407 Nasal Mask



Presented With Compliments From
www.DirectHomeMedical.com

The **FlexiFit™ 407** Nasal Mask, developed using our anthropometric database with thousands of facial contours, provides a better fit and superior comfort while maintaining the ease of use features common to the **FlexiFit™** Series.


The **FlexiFit™** Technology has been adapted from the **FlexiFit™ 405** to provide a greater range of movement at the nasal bridge, simplifying the fitting process while improving the seal and comfort. The unique **Glider™** Strap remains for freedom of movement while maintaining a seal. The new **Stretchgear™** Headgear maximizes comfort and minimizes facial pressure during head movement.

The **FlexiFit™ 407** is designed and contoured especially to meet the goal of one size to fit all, and comes completely assembled for a complete and simple Out-of-Box Solution.

Product Code	Description	Quantity
HC407	FlexiFit™ 407 Nasal Mask includes: <ul style="list-style-type: none">• Silicone Seal (one size)• Headgear• Glider™ Strap	1 1 1


BREATH-O-PRENE™ is a registered trade mark of AccuMED Technologies, Inc

FlexiFit™ 407 Mask Spare Parts



Elbow Kit

Product Code	Description	Quantity
400HC202	Elbow Non-Diffuser (Holes) plus Swivel	1 of Each
400HC200	Elbow Non-Diffuser (Holes) plus Swivel	10 of Each
900HC441	Elbow Plus Diffusers material, Swivel	1 of Each



Foam Cushion

Product Code	Description	Quantity
400HC002	Foam FlexiFit™ 407	Each




Foam Cushion & Silicone Seal Combo

Product Code	Description	Quantity
400HC501	Single Foam & Silicone Seal FlexiFit™ 407	Each




Swivel Connector

Product Code	Description	Quantity
900HC438	Swivel Connector Female	10 Pack



Oxygen/Pressure Port Connector

Product Code	Description	Quantity
900HC452	Oxygen/Pressure Port Connector	10 Pack




Forehead Pads

Product Code	Description	Quantity
900HC423	Forehead Pads - fits FlexiFit™ 405, 407, 431	2 Pack




Glider

Product Code	Description	Quantity
900HC418	Glider™ Strap - fits FlexiFit™ 405 and 407	10 Pack



Mask Headgear

Product Code	Description	Quantity
400HC301	Stretchgear™ Headgear FlexiFit™ 407	Each



Mask Kit

Product Code	Description	Quantity
400HC502	FlexiFit™ 407 Mask – No Headgear	Each