


Direct Nasal Mask

INFINITY™ 481 Direct Nasal Mask



Presented With Compliments From
www.DirectHomeMedical.com

The **Infinity™ 481** Direct Nasal Mask is smaller and lighter than all other masks to offer the patient total sleeping freedom.

The **Infinity™ 481** differs from other over-the-nose nasal masks by delivering pressure directly into the nasal passages via an anatomically shaped Silicone Seal direct to each nostril. A maintenance free vent disperses the air quietly and efficiently, while a ball and socket styled joint allows infinite movement for comfort.

This freedom of movement, combined with minimalist **Microgear™** Headgear provides for unprecedented comfort and ease of use. A conventional headgear option is also included for additional stability.

Product Code	Description	Quantity
HC481	Infinity™ 481 Direct Nasal Mask includes: <ul style="list-style-type: none">• Mask Base and Flexitube Assembly• Microgear™ Headgear• Headgear• Silicone Seal (XS, S, M, L sizes)• Sizing Guide• Storage Holder	1 1 1 1 1 1

Infinity™ 481 Mask Spare Parts

Silicone Seal




Product Code	Description	Quantity
400HC107	Silicone Seal Infinity™ 481 XS, S, M, L Plus Sizing Guide	1 of Each
400HC108	Silicone Seal Infinity™ 481 XS	Each
400HC109	Silicone Seal Infinity™ 481 S	Each
400HC110	Silicone Seal Infinity™ 481 M	Each
400HC111	Silicone Seal Infinity™ 481 L	Each



Headgear

Product Code	Description	Quantity
400HC303	Infinity™ 481 Microgear™ Headgear	Each
400HC308	Infinity™ 481 Headgear BREATH-O-PRENE™ Plus Microgear™ Headgear	1 of Each



Sizing Guide

Product Code	Description	Quantity
400HC511	Sizing Guide	10 Pack



Mask Base

Product Code	Description	Quantity
400HC505	Infinity™ 481 Mask Base (Base, Tube, Connectors)	1 of Each