


 Sleepnet®  
Fine, handcrafted sleep masks

**Mojo™**

**Full Face Mask  
Size Gauge**




Hold size gauge up to patient's face and measure from indent on nose bridge to indent under lips. Relax jaw, with mouth slack and slightly open (about 1 1/2.5cm).



Presented With Compliments From  
[www.DirectHomeMedical.com](http://www.DirectHomeMedical.com)

90102 Rev. B



Print out template at 100%.

Make sure not to choose "shrink to fit" when printing.

To ensure your template is the correct size, measure the "1 Inch" square above.