

Oral Mask

ORACLE™ 452 Oral Mask



Presented With Compliments From
www.DirectHomeMedical.com

The Fisher & Paykel Healthcare Oracle™ was the first ever developed oral mask, a milestone in interface design.

The **Oracle™ 452** delivers pressure through the mouth rather than the nose offering the patient yet another choice of delivery. Clinically proven to function as effectively as a nasal mask, the **Oracle™ 452** is a perfect solution for mouth breathers and patients with chronic nasal obstruction who find therapy difficult with a nasal mask.

The new generation **Oracle™ 452** has exciting refinements that improve the patient experience, for ease of use and ease of seal. A new symmetrical design means the mask is always correctly fitted while the infinitely adjustable rotating dial allows for a personalized fit to provide greater comfort and simplicity.

The **Oracle™ 452** Interface is a new, unique oral mask for the delivery of CPAP or bi-level therapy.

Product Code	Description	Quantity
HC452	Oracle™ 452 Oral Mask includes: <ul style="list-style-type: none">• Flexitube Assembly• SoftSeal (small)• SoftSeal (large)• Headgear with Easy-Release Cord• Storage Case	1 1 1 1 1

Oracle™ 452 Mask Spare Parts



Diffusers

Product Code	Description	Quantity
900HC439	Diffuser Kit also fits FlexiFit™ 405	10 Pack



Diffuser Caps

Product Code	Description	Quantity
900HC440	Diffuser Caps also fits FlexiFit™ 405	10 Pack




Elbow Kit

Product Code	Description	Quantity
900HC441	Elbow Kit also fits FlexiFit™ Includes material, cap and swivel	1 of Each



Silicone Seal

Product Code	Description	Quantity
400HC105	Oracle™ 452 SoftSeal Small	Each
400HC106	Oracle™ 452 SoftSeal Large	Each



Oxygen/Pressure Port Connector

Product Code	Description	Quantity
900HC452	Oxygen/Pressure Port Connector	10 Pack



Flexitube Kit

Product Code	Description	Quantity
400HC209	Oracle™ 452 Flexitube Kit, Elbow & Vented Non Re-Breathing Valve	Each
400HC205	Oracle™ 452 Flexitube Kit Tube, Elbow & Non Re-Breathing Valve (Australia Only)	Each

Oracle™ 452 Mask Spare Parts

Mask Headgear



Product Code	Description	Quantity
400HC300	Oracle™ 452 Headgear Accessory plus clips	Each
400HC215	Oracle™ 452 clips only	2 Pack

Mask Kit



Product Code	Description	Quantity
400HC506	Oracle™ 452 - No Headgear	Each

Nasal Plugs



Product Code	Description	Quantity
400HC204	Oracle™ 452 Nose Plugs S,M,L	1 of Each

Non Re-Breathing Valve



Product Code	Description	Quantity
400HC206	Vented Non Re-Breathing Valve	Each
400HC208	Non Re-Breathing Valve (Australia Only)	Each