




Ultra Mirage™ II


NASAL MASK

Mirage Activa™


NASAL MASK




↑
STANDARD




↑
LARGE



↑
SHALLOW



↑
SHALLOW WIDE
(Ultra Mirage II only)



Presented With Compliments From
www.DirectHomeMedical.com

How to use the fitting template:

1. Cut out the template along the dotted line
2. Cut out nasal holes
(hint: fold each nasal hole in half to cut)
3. Follow directions on template to determine appropriate cushion size

Actual Size: 5 1/2" (Width) x 7" (Height)

Ruler is for scaling purposes.

RESMED