




The Essential  
"4th element"

fits within any CPAP machine travel case



'Folding, yet all parts strung gathered in order, it weighs only 5.5 oz (156g) and is a very small 2.1" diameter x 7" long.

Dual Use: Home & Travel

The new improved CPAP Hose Lift v5.5 is simply the lightest, most compact home and travel, Dual Use hose management device available.


Compared to Previous v5

12% Lighter

16% More Compact


14% Larger Base

New: Optional Height / Reach



**NEW**  
with v5.5

**Optional height and reach – just flip the bend, and then the v5.5 will reach out a bit farther and stand a bit shorter. Start by 'unknotting' bungee at arrow:**



CE mark: CA011750; HTS#: 9019.20.00;  
NPI#: 1467684712; Insurance Code#: 1399; Model: CPAPHLv5.5



The first and (only) true *Original*, with a Curly Cord "roller" on the suspending loop, eases movement of the cpap hose within.



A deep socket mast connection, a mast to base [stretched] cord, "tighten-n-lock:" mast to base - provides an even more stable platform, in hose suspension support.



Includes a Hook for your mask to store, out of your way.