

Review Safety Checklist Before Using Your SoClean 3

The following steps **must** be completed to setup and activate your device:

- Sleep equipment accessories are washed in sink before first use (page 20)
- SC3 is powered on (page 12)
- SoClean Adapter is installed on sleep equipment machine (page 14)
- Ozone tube is connected to sleep equipment machine (page 16)
- Sleep equipment hose is still connected to sleep equipment machine and mask
- Filter is installed (page 17)
- Device Unlock/Activation

Press the Control Panel buttons in the following order to activate your SoClean 3:



For more information on SoClean 3, visit [DirectHomeMedical.com](https://www.DirectHomeMedical.com)
For additional assistance, call **888-505-0212**