



Designed with compliance in mind

InnoSpire Deluxe compressor nebulizer system

For use with aerosol medications, InnoSpire Deluxe incorporates built-in conveniences that encourage adherence to aerosol medication treatment plans. As part of the InnoSpire family of compressor nebulizer systems, Deluxe combines a durable compressor in a modern design that provides delivery of consistent aerosol medication treatments.

InnoSpire Deluxe is designed to include a convenient storage compartment to hold extra medicine or nebulizer supplies. The attached power cord tucks neatly into the storage area and a sturdy handle makes it easy to carry. It is covered by a five-year warranty and is compliant with the international basic safety and essential performance standard IEC 60601-1 third edition and its amendments.



Presented With Compliments From
www.DirectHomeMedical.com

PHILIPS
RESPIRONICS

InnoSpire Deluxe



Key features

- Robust compressor suitable for continuous use
- Integral storage for nebulizer accessories
- Lid design incorporates 180° hinge and secure latch and can be positioned flat for additional surface space
- Neb holder provides an extra 'hand' for stability when filling the cup or to park the handset if taking a rest during treatment



Venturi system

provides a faster delivery rate resulting in more effective treatment times than other conventional nebulizers

Innovative design consistently produces high-quality aerosol

InnoSpire Deluxe partners with SideStream nebulizers to provide fast and efficient drug delivery. SideStream nebulizers feature an active venturi system that provides air flow in addition to that supplied by the compressor, resulting in improved drug delivery and more effective treatment times.



SideStream Plus

Product specifications	
Size	380mm x 177mm x 110mm (13" D x 6.7" W x 4.3" H)
Weight	2.0 kg (4.4 lbs)
Volume level	58±3dBA
Electricity supply	120V/60Hz, 1.7 Amp
Thermal overload protector	Resettable thermal fuse; functioning temperature 130C°
Maximum pressure	268 kPA (35 psi)
Flow rate	6 LPM @ 10 psi
Max flow rate	10.78 LPM
Particle size distribution	0.5 to 5 microns

Ordering information			
Product description	Reorder number	HCPCS code	NDC number
InnoSpire Deluxe compressor with SideStream disposable nebulizer (6/cs)	1110063	E0570	08373-0063-00
Filter (4 pack)	1102088	A7013	08373-2088-00
SideStream disposable nebulizer (50/cs)	HS800	A7003	08373-0800-00
SideStream reusable nebulizer (10/cs)	HS860	A7005	08373-0860-00
SideStream Plus nebulizer (10/cs)	HS870	A7005	08373-0870-00
SideStream pediatric mask (50/cs)	1025529	A7015	08373-5529-00
SideStream adult mask (50/cs)	1025531	A7015	08373-5531-00
SideStream Plus adult mask (50/cs)	1052521	A7015	08373-2521-00
Sami the Seal child-friendly mask (50/cs)	HS880	A7015	08373-3461-00

- Reusable SideStream is recommended for up to six months of use
- Disposable SideStream is recommended for up to two weeks of use

Please visit www.philips.us/innospire-deluxe for more information

Respironics, InnoSpire Deluxe and SideStream Plus are trademarks of Respironics, Inc. and its affiliates.

© 2014 Koninklijke Philips N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Geyer JH 12/20/13 MCI 4105959 PN 1114498

Philips Healthcare is part of Royal Philips

Philips Respironics International Headquarters
+33 1 47 28 30 82
Philips Respironics Asia Pacific
+65 6882 5282
Philips Respironics Australia
+61 (2) 9947 0440
Philips Respironics United Kingdom
+44 800 1300 845

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)
www.philips.com/respironics

