



DESIGNED WITH **HEAVY USE**

PATIENTS IN MIND



Presented With Compliments From
www.DirectHomeMedical.com



AEROSOL DELIVERY SYSTEM



The Vios® PRO Aerosol Delivery System was designed with heavy use patients in mind. Available with both PARI LC PLUS® and LC® Sprint nebulizers, the powerful 1.6 bar compressor withstands repeated use and provides fast treatments.

www.pari.com/viospro

TIME TO GO **PRO!**



Configurations

	Part #	HCPCS #
Vios PRO, LC Sprint	312F35-LCS	E0570, A7005
Vios PRO, LC PLUS	312F83-LC+	E0570, A7005

Technical Specifications

Operating Pressure ³	23.2 PSI (1.6 bar)
Liter Flow ³	5.1 LPM
Weight	4 lbs
Dimensions (L, W, H)	5.75 x 7.5 x 5.5 in
Compressor Warranty	5 years

Performance with LC[®] Sprint and LC PLUS[®]

LC Sprint¹

Total Output Rate:	590 mg/min
MMD (Mass Median Diameter)	2.9 µm
RF% (Mass Percentage <5 µm)	75%

LC PLUS²

Total Output Rate:	440 mg/min
MMD (Mass Median Diameter)	3.7 µm
RF% (Mass Percentage <5 µm)	67%

¹Measurement with the Malvern Mastersizer X at 23°C, 50% relative humidity, nebulized medium: 0.9% NaCl (5ml), inspiratory flow 20l/min.

²Measured with Malvern MasterSizer X at 23° C, 50% relative humidity, 0.9% NaCl solution, inspiratory flow 20 liters/minute, continuous nebulization, fill volume 2.5 ml.

³When used with a PARI LC[®] reusable nebulizer.

PATIENTS **BENEFIT** FROM:



Fast treatments allow patients to spend less time on aerosol therapy helping to improve their quality of life



Available with the clinically proven PARI LC PLUS[®] or LC[®] Sprint nebulizers



Powerful 1.6 bar compressor withstands repeated use for heavy use patients

Accessories



PARI LC SPRINT
023F35-VP



Replacement Air Filters
085F0012P2



PARI LC PLUS
022F81-VP



Vios[®] PRO Tubing
041F4594

Note: Only use Vios[®] PRO tubing with the Vios[®] PRO Aerosol Delivery System



PARI Respiratory Equipment, Inc.
Tel: 1-800-FAST-NEB (327-8632)
product.info@pari.com
www.pari.com | Midlothian, VA 23112