



**ResMed**

ResMed **Air Solutions**

great sleep



DirectHome  
MEDICAL

Presented With Compliments From  
[www.DirectHomeMedical.com](http://www.DirectHomeMedical.com)



a new beginning  
in great sleep

[ResMed.com/AirSolutions](http://ResMed.com/AirSolutions)

Introducing a complete sleep apnea therapy solution specifically for women. ResMed Air Solutions combines the first therapy device designed to provide female-specific therapy with ResMed's for Her masks, cutting-edge technology and tailored support to provide everything you need for a great night's sleep.





## Out-of-the-box comfort

The AirSense™ 10 AutoSet™ for Her has been designed to make therapy easy. With no complicated settings or menus to navigate, simply press Start and you're ready to go. Each device comes with the HumidAir™ heated humidifier built in so you can enjoy the benefits of humidification in one easy-to-use system. And when used with the optional ClimateLineAir™ heated tube, your device has been designed to automatically deliver the optimal temperature and humidity for the ultimate comfort experience.

## Unique therapy for women

By responding to female-specific characteristics of sleep apnea, ResMed's new AutoSet for Her algorithm has been designed to meet your individual needs, providing a tailored approach to sleep therapy. And with the new ResMed AirFit™ for Her mask series, it's easier than ever to find the right fit for you.

Choose from the AirFit P10 for Her nasal pillows mask, the AirFit N10 for Her nasal mask or the AirFit F10 for Her full face mask for a comfortable, stylish solution designed to make therapy easy—no matter which one you choose.



## A comfortable start

The new AutoRamp™ feature with sleep onset detection helps make therapy comfortable from the moment you turn on your device. AutoRamp works by delivering a low pressure to help you fall asleep with ease. Once it detects you're asleep, it comfortably increases the pressure to ensure your prescribed level is delivered the moment you need it.


## A restful night

Every device has an ambient light sensor that adjusts the brightness of the screen based on the amount of light in the room. So if you get up during the night and need to pause your therapy, the screen will stay dim so it doesn't disrupt your partner. The screen's back light also turns off automatically to help you fall back to sleep without disturbance. And with ResMed's Easy-Breathe motor built into each device, you and your partner can enjoy whisper-quiet sleep, every night.

Find out more at  
[ResMed.com/AirSolutions](https://www.ResMed.com/AirSolutions)

## Personalized support

Like your own sleep coach, ResMed's new online program, myAir,<sup>™</sup> provides support during every stage of your treatment journey. By giving you access to your therapy data, it lets you continually track your progress and see how well you're doing. Emails or text messages can also provide helpful education and encouragement to support you along the way. As you mark your successes and celebrate your milestones with awards and badges, you'll stay motivated to keep going.





Distributed by ResMed Corp, 9001 Spectrum Center Boulevard, San Diego, CA 92123 USA +1 858 836 5000 or 1 800 424 0737 (toll free). See ResMed.com for other ResMed locations worldwide. AirFit, AirSense, AutoRamp, HumidAir and myAir are trademarks of ResMed Ltd. AutoSet and ClimateLine are trademarks of ResMed Ltd and are registered in the U.S. Patent and Trademark Office. Specifications may change without notice. ©2014 ResMed. 1017633/1 2014-06

**R<sub>x</sub>** Only

[ResMed.com/AirSolutions](http://ResMed.com/AirSolutions)