



Presented With Compliments From
www.DirectHomeMedical.com



DeVilbiss IntelliPAP® Standard CPAP System

You talked. We listened.
Everyone gets what they need.

The new **DeVilbiss IntelliPAP Standard CPAP System** has been designed with patients and providers in mind, incorporating many of the recommendations solicited through research to **optimize patient adherence**. IntelliPAP's combination of **comfort, education and adherence tracking with SmartCode™** helps to ensure patient compliance – all in a **highly efficient platform** that makes great business sense to you.

From its **reliable electronics and water resistant housing** to its durable Lexan® humidifier chamber and stainless steel heat transfer plate, IntelliPAP is **built for years of dependable performance and patient satisfaction**. Just what you would expect from DeVilbiss – the widely recognized leader in respiratory technology.

DeVilbiss IntelliPAP
Standard CPAP Unit
Model DV51D-HH

Performance Features

- Quietest **CPAP Unit on the market** – ONLY 26 dBA!*
- Integrated Heated Humidification
- **Automatic Altitude Adjustment** up to 9,000 feet
- Exclusive humidifier chamber design **maximizes humidification and minimizes leaks**
- Exceptionally small footprint, **smallest integrated unit** among top-selling CPAP
- Delay pressure ramp
- Best-in-industry warranty



Try our Serenity®, FlexSet® or FlexAire® masks – DeVilbiss offers a lightweight, comfortable and latex-free line of masks. Visit www.DeVilbissHealthcare.com for more information.

Design Features

- Clean, sophisticated design
- Whisper-quiet operation – 26 dBA*
- Exceptionally small footprint fits easily into the home environment and carry-on luggage
- Discreet carry bag included that can be easily expanded to hold additional accessories
- Integrated heated humidification system easily docks into the bottom of the unit without requiring a separate hose or power cord
- Humidifier reservoir is easy to remove and replace and includes built-in, overflow protection
- Water ingress protection prevents humidifier water from damaging the CPAP – great for travelers
- Patient air tubing extends from back of unit to reduce clutter and enhance visual appeal
- Accommodates various tubing lengths (6 and 10 feet)
- Large, slip-resistant pads on bottom of unit to minimize movement on night table during the night
- Eight-foot power cord length provides greater positioning flexibility – no bulky power brick on cord
- Easy and intuitive set up – three, large buttons control nearly all functions for simplified patient operation
- Large, backlit LCD and easy-to-follow patient menu
- Auto on/off automatically starts the flow of air after breathing into the mask and stops the flow of air when mask is removed
- Delay pressure ramp eases the machine to the prescribed pressure – allowing the user to fall asleep at a lower pressure (0-45 minutes in 5 minute increments and .5 increments for pressure)
- SmartCode™ compliance tracking encodes usage data that can be communicated to a clinician and then translated
- Reminders alert to the need to replace or check parts
- Extra air inlet and fine particle filter included
- Industry leading, 3-year warranty

Patient Education and Marketing Support

- User video includes an introduction to CPAP therapy with tips for acclimating to, and optimizing, therapy, as well as how to care for and use the IntelliPAP system
- Patient instruction guide demonstrates how easy it is to use IntelliPAP, while providing a comprehensive overview of sleep apnea, its treatment and tips for acclimating to therapy



DeVilbiss IntelliPAP Standard CPAP System includes attractive travel bag

Technical Specifications

* Tested per ISO 17510

| | | | |
|--------------------------------------|-------------------------------|--|---|
| Product Dimensions | 6.4" H x 6.5" W x 8.4" D | Automatic Leak Compensation | Yes |
| Product Weight | 4.45 lbs w/humidifier and PAP | Auto On/Off | Yes (can be enabled or disabled) |
| Pressure Range | 3-20 cm H2O (.5 increments) | Motor Hour Meter | Yes |
| Universal Power Supply | Yes | Therapy Compliance Meter | Yes |
| AC Cord Length | 8 feet | Warranty | 3-year warranty on IntelliPAP 2-year warranty on humidification system |
| Automatic Altitude Adjustment | Yes | DeVilbiss SmartLink™ Compatible | Yes |
| DC Operable | Yes | SmartCode™ Feature | Yes |
| ABOB Analog Output for PSG | Yes (pressure only) | | |



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com