Sleep Reporting Made Easy

Put your Compliance Data in the Palm of Your Hand.



- For Android and Apple devices
- Compatible with all Luna® models
- As simple as taking a picture free of charge

The first of its kind

~ 3B Luna® QR app lets you take your compliance data with you.



REPORTING IN THREE EASY STEPS

Step 1:

Search for 3B[™] Luna[®] QR in the Google Play or App Store and download the FREE app.

Step 2:

Once installed, link the app to your Luna PAP device.

Step 3:

Scan the iCodeQR or iCodeQR+ that is generated on your PAP device. In seconds your summary compliance data is saved in your 3B™ Luna® QR app and is accessible anywhere you need it.

COMPLIMENTARY SUPPORT

Contact 3B Medical directly by visiting us at www.3Bproducts.com and clicking on "Support".

QUICK REPORT FEATURE IN ICODECONNECT™

Step 1:

Go to www.icodeconnect.com

Step 2:

Click on the "Quick Report" button.

Step 3:

Enter required field information; name and report type



Add: 203 Avenue A NW, Suite 300 | Winter Haven FL 33881 Tel: 863-226-6285 Fax: 863-226-6284