



At home or on-the-go,
Transcend provides the
sleep apnea therapy you need.

Featuring its small, light and innovative
design, Transcend is the gold standard
in portable PAP innovation. With a full
line of PAP products and accessories
unlike any other CPAP on the market,
you can customize your Transcend
device to fit your active lifestyle. Sleep
apnea therapy is no longer limited to
the night stand.



SLEEP APNEA THERAPY

Customized therapy
is as easy as 1-2-3

Transcend's patient-centric product line
features accessory options unlike any other
sleep apnea therapy system on the market
to give you therapy that fits your lifestyle.

1 Therapy

- Transcend Auto™
with expiratory relief
- Transcend EZEX™ CPAP
with expiratory relief
- Transcend® CPAP
- TranSync™ wireless therapy monitoring



2 Humidification

- Transcend Heated Humidifier™



3 Power Options

- P8™ Battery *
- P4™ Battery *
- Transcend Portable
Solar Battery Charger™
- DC Mobile Power Adaptor



* Not for use with Transcend Heated Humidifier.

SLEEP APNEA THERAPY



Sleep apnea therapy
for the active lifestyle.



Presented With Compliments From
www.DirectHomeMedical.com



www.MyTranscend.com © 2013 Somnetics International, Inc.
Protected by U.S. and foreign patents issued and pending. Transcend® is a registered
trademark of Somnetics International, Inc. - PN: 400130 Rev. B

Transcend is changing the way sleep apnea therapy is delivered.™

Identical in size, each Transcend device is designed to fit your active lifestyle.



Transcend AUTO™

Transcend Auto is a fully-featured APAP that continually monitors your breathing to adapt to your changing therapy needs. Transcend EZEX technology provides pressure relief on exhalation for more breathing comfort.



Transcend EZEX™

Transcend EZEX is a fully featured fixed-pressure CPAP with EZEX exhalation pressure relief technology for added comfort.



Transcend®

Transcend is a fully featured fixed-pressure CPAP. This original CPAP is the pioneer in portable sleep apnea therapy. At home or on the road, Transcend is designed to enhance the lifestyles of today's active CPAP users.

Small size, big features for therapy anywhere.

Transcend's fully-featured therapy devices are small, lightweight and portable. They are perfect for packing into a suitcase or taking up minimal space on the night stand.



Items shown proportionate.

	Transcend AUTO™	Transcend EZEX™	Transcend®
Auto-titrating	✓		
Fixed pressure	✓	✓	✓
Operating pressure (4-20 cmH2O)	✓	✓	✓
Transcend EZEX™ pressure relief (settings 1,2,3)	✓	✓	
Ramp (0-45 min), adjustable	✓	✓	✓
Weights less than one pound (426 g)	✓	✓	✓
Automatic altitude adjustment	✓	✓	✓
AHI & leak detect	✓	✓	✓
Up to 13 months compliance data storage	✓	✓	✓
TranSync™ Wireless Compliance capability	✓	✓	
Compatible with Transcend batteries	✓	✓	✓
Compatible with Transcend Heated Humidifier™	✓	✓	✓

For full specifications, please visit www.myTranscend.com

Contact your medical supplier for more information.