

Economies of scale



Presented With Compliments From
www.DirectHomeMedical.com

PROBASICS™

by PMI

Zzz-PAP™



Introducing the new **Zzz-PAP™** CONTINUOUS POSITIVE AIRWAY PRESSURE UNIT

The Probasics Zzz-PAP™ by PMI delivers performance, simplicity, style, and value in a very small package. The Zzz-PAP offers easy patient setup with a 4-18 cmH₂O pressure range, a comfortable 0-45 minute ramp feature and manual altitude adjustment. PMI also offers the optional Probasics Zzz-PAP Heated Humidifier that quickly and easily integrates to the Zzz-PAP unit providing excellent relief from the symptoms of dryness and congestion associated with CPAP therapy. Weighing less than two pounds and coming standard with an auto-adjusting 100 Volt/240 Volt Universal Power Supply, the compact Zzz-PAP is not only an excellent choice for your bedside but it's perfect for travel as well. Best of all, the Probasics Zzz-PAP is backed by a two-year manufacturer's warranty which we all know is the most necessary feature for providing both the patient and the dealer with peace of mind.



7501

Zzz-PAP™ INTEGRATED HEATED HUMIDIFIER

Designed for use with the Probasics Zzz-PAP™ Mini-CPAP unit, PMI offers an Integrated Heated Humidifier system to relieve nasal irritation and dryness which result from constant airflow by supplying the required humidity to the upper airway. With a reliable and safe design which features overheat protection, the Probasics Zzz-PAP Heated Humidifier can drastically improve the patient's comfort level and improve patient compliance.



7600



7500

7500	Zzz-PAP CPAP Unit with Interated Heated Humidifier
7501	Zzz-PAP CPAP Unit alone
7601	Zzz-PAP Interated Heated Humidifier alone

SPECIFICATIONS



Continuous Positive Airway Pressure

Easy-to-carry design

It is compact and easy to carry with its small size and light weight...Perfect for travel.

LCD display with back light and easy-to-use interface

The LCD display is easy to read, and the intuitive user interface is very simple to use.

Easy Pressure level adjustment

Operating pressure, ranging from 4 to 18 cm H₂O, can be easily and securely set to comply with the physician's prescription while in provider mode.

Easy Ramp time setting and adjustment

Deluxe ramp capability allows a 0 - 45 minute adjustable time period to help user fall to sleep easier.

Compliance meter

The number of patient use hours can be easily accessed for reimbursement or therapeutic recording purposes.

Extremely quiet operation

Performs at just 30 dBA at 10.0 cmH₂O.

Pressure Range	4.0 - 18.0 cmH ₂ O (0.5 cmH ₂ O increment)
Pressure Accuracy	+2 -1 cmH ₂ O
Ramp Time	0 - 45 min (5 minute increments)
Altitude Compensation	Level 1 ~ 8 for 0 ~ 8000 ft. (manual setting)
Device Operating Meter	9999.9 hr
Heated Humidifier	Yes
Sound Level	30 dBA @ 10 cmH ₂ O, one meter distance
Dimension (W x D x H)	14.5 x 13.0 x 10.0 cm 5.7" x 5.1" x 3.9"
Weight	800 g (estimation)
Power Supply	100V-240V, 50/60Hz, Universal Power Supply
Operating Temperature	5°C to 35°C / 41°F to 95°F
Environment Humidity	15% to 95% RH



Integrated Heated Humidifier

- Integrates with Probasics Zzz-PAP mini CPAP unit
- Easy to use rotary knob for heating level adjustment
- Dishwasher safe water chamber has an easy to remove base plate for cleaning
- Safe & reliable with overheat protection
- Indicators clearly show when Power is ON and Heater is ON

Power Supply	120VAC, 50/60Hz 230VAC, 50/60Hz
Water Capacity	450 ml
Heater Settings	1 - 6 (Heater plate temperature 86°F to 140°F / 30°C to 60°C)
Warm-Up Time	< 30 minutes
Pressure Drop	0.2 cmH ₂ O @ 60LPM
Humidity Output	10 - 40 mgH ₂ O/L
Indicator	Power On indicator, Heater On indicator
Dimension (W x D x H)	Heater: 6.8" x 10" x 1.6" (17.2 x 25.5 x 4 cm) Water Chamber: 4.5" x 5.9" x 4.2" (11.5 x 14.9 x 10.6 cm)
Weight	Heater: 380 g (0.84 lb) Water Chamber: 240 g (0.53 lb)

All specifications are subject to change without notice.

0607



Probasics™

15 South Main Street

Second Floor, Unit F

Marlboro, NJ 07746

Tel. 800.438.0503

Fax 732.683.1945

www.pmiusa.biz