

ResMed





Presented With Compliments From
www.DirectHomeMedical.com

Use this fitting template as a guide to find the correct size for all AirFit N10, AirFit N10 For Her, Swift FX Nano and Swift FX Nano For Her Masks.

To use the fitting template:

1. Cut out the template along the dotted line
2. Follow directions on template to determine the recommended cushion size.
3. Please ensure the template is printed to scale using the ruler at the bottom of this sheet.



AirFit and Swift are trademarks of ResMed Ltd. Swift is registered in the U.S. Patent and Trademark Office.
 © 2014 ResMed Ltd. All rights reserved. Specifications may change without notice. Product availability may vary across regions.
Rx Only  0123

Ruler is for scaling purposes.

