



## Upgrade to blue

**ComfortGel Blue Full** is an excellent advancement for patients who need a full-face mask.

Comfort and design make the ComfortGel Blue Full a sensible choice for more of your patients.

A comfortable forehead cushion and enhanced StabilitySelector combine with our popular blue gel replaceable cushion to create a high level of comfort and stability. The unique sure seal technology (SST) provides a membrane-thin layer over the gel cushion to closely follow facial contours of both men and women, allowing more

wearers to quickly achieve the comfortable seal they need. With ComfortGel Blue Full, your full-face mask patients finally have a simple and sensible way to enjoy the restful sleep they need. It's one more way we're helping you fit 100% of your patients, 100% of the time.




Presented With Compliments From  
[www.DirectHomeMedical.com](http://www.DirectHomeMedical.com)

**PHILIPS**

**RESPIRONICS**

sense **and** simplicity



- Premium forehead pad enhances comfort and fit.
- Premium headgear provides patient comfort and mask stability.
- StabilitySelector moves freely and is easy to adjust.
- O<sub>2</sub> pick-off port for pressure measurement.
- Angled exhalation micro ports make operation quieter and redirect air away from a bed partner.
- Ball-and-socket Quick Clips mean freer movement and easy mask removal.

- Integrated retaining ring tabs simplify removal.
- Replaceable gel cushion creates comfort and seal.
- Silicone flap with sure seal technology (SST) optimizes seal.

Item	Small	Medium	Large	Extra large	U.S. HCPCS code
<b>ComfortGel Blue Full Mask</b>					
With headgear	1081800	1081801	1081802	1081803	A7030, A7035
Mask only	1081810	1081811	1081812	1081813	A7030
DuoPack: mask and headgear plus 2 same-size cushions	1081805	1081806	1081807	1081808	A7030, A7035, A7031
DuoPack: mask only plus 2 same-size cushions	1081815	1081816	1081817	1081818	A7030, A7031
<b>ComfortGel Blue Full Replacement cushions/flaps</b>					
Cushion and SST flap	1081895	1081896	1081897	1081898	A7031
Cushion only	1081890	1081891	1081892	1081893	A7031
<b>Parts and accessories</b>	<b>Part number</b>				<b>U.S. HCPCS code</b>
ComfortGel Blue Full headgear	1040138				A7035
RS Premium headgear (smaller)	1059236				A7035
Forehead pad	1040114				N/A
Entrainment valve swivel	1081899				N/A
ComfortGel Blue Full Quick Clips (2 per pack)	1040139				N/A
Port cap (5 per pack)	1061825				N/A

Respironics, ComfortGel, DuoPack, Quick Clip, SST, and StabilitySelector are trademarks of Koninklijke Philips Electronics N.V. and its affiliates.

Please visit [www.philips.com/respironics](http://www.philips.com/respironics)



© 2012 Koninklijke Philips Electronics N.V. All rights are reserved. Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication. CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Hoehc JJ 2/08/12 MCI 4104739 PN 1099038

Philips Healthcare is part of Royal Philips Electronics

Philips Respironics  
1010 Murry Ridge Lane  
Murrysville, PA 15668

Philips Respironics International  
Headquarters  
+33 1 47 28 30 82

Customer Service  
+1 724 387 4000  
800 345 6443 (toll free)  
[www.philips.com/respironics](http://www.philips.com/respironics)