

## Tried and true

**ComfortGel Blue** – the next generation of comfort and technology, delivering everything you've come to expect.

An even higher level of patient comfort and an enhanced overall experience make ComfortGel Blue a mask choice more of your patients will embrace.

Staying true to our time-proven blue gel masks, we've retained the globally popular comfort and flexibility features you expect, while enhancing patient comfort even more. We've incorporated an improved forehead pad designed to help reduce pressure points,

and developed a lower profile exhalation port with an integrated swivel that quietly directs air flow up and away from a bed partner. We've also adopted a new gel formulation. All designed with an enhanced patient experience – and better compliance – in mind. And, using ComfortGel Blue with Philips Respironics System One Resistance Control will deliver optimum PAP therapy and therapy comfort. ComfortGel Blue and our complete line of masks go farther than ever to help you fit 100% of your patients, 100% of the time.

**PHILIPS**


**RESPIRONICS**

One-size-fits-most premium headgear with EZ Peel tabs adjusts easily in four places for a secure, custom fit.

Silicone Comfort Flap improves seal and can be removed for easy cleaning.

New generation gel technology cushion is soft and conforms to patient's face for added comfort.

Ball-and-socket headgear attachments enable free movement and easy mask removal.




Premium forehead pad improves comfort and fit.

StabilitySelector raises and lowers for proper mask angle, seal, and fit.

Angled exhalation micro ports assure quiet operation and redirect air away from a bed partner.

Two 360° swivels allow tubing to move more freely.

Product name and part number are included on faceplate for easy identification by you and patients.



Sizes and System One Resistance Control settings are printed on bottom of cushion.

Ordering information					
ComfortGel Blue mask	Petite	Small	Medium	Large	U.S. HCPCS
Mask with headgear	1070040	1070039	1070038	1070037	A7034, A7035
Mask only	1070044	1070043	1070042	1070041	A7034
DuoPack (includes headgear and two cushions of the same size)	1070052	1070051	1070050	1070049	A7034, A7035, A7032
FitPack: includes headgear and small and medium cushions	N/A	1070053		N/A	A7034, A7035
Replacement cushions	Petite	Small	Medium	Large	U.S. HCPCS
ComfortGel Blue cushion and flap	1070108	1070107	1070106	1070105	A7032
ComfortGel Blue cushion	1070104	1070103	1070102	1070101	A7032
Replacement parts			Part number	U.S. HCPCS	
Premium forehead pad			1040114	N/A	
Premium headgear			1033678	A7035	
Ball-and-socket Quick Clips, 2 per pack			1009723	N/A	
ComfortGel Blue sizing gauge			1012040	N/A	

Respironics, ComfortGel, Comfort Flap, DuoPack, FitPack, Quick Clip, StabilitySelector, and System One are trademarks of Respironics, Inc. and its affiliates, all rights reserved.

Please visit [www.philips.com/comfortgelblue](http://www.philips.com/comfortgelblue)



© 2010 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Hoeh JJ 4/21/10 MCI 4103170 PN 1071163

Philips Healthcare is part of Royal Philips Electronics

Respironics Asia Pacific  
+65 6298 1088  
Respironics Australia  
+61 (2) 9666 4444  
Respironics Europe, Middle East, Africa  
+33 1 47 52 30 00  
Respironics United Kingdom  
+44 800 1300 845

Philips Respironics  
1010 Murry Ridge Lane  
Murrysville, PA 15668

Customer Service  
+1 724 387 4000  
800 345 6443 (toll free, US only)  
[www.philips.com/respironics](http://www.philips.com/respironics)