



Presented With Compliments From  
[www.DirectHomeMedical.com](http://www.DirectHomeMedical.com)

## DeVilbiss Serenity<sup>®</sup> Nasal CPAP Mask

### Simple and Intuitive Design

The simple and intuitive design of this economy nasal CPAP mask makes it the perfect choice. The Serenity mask features easy to size headgear and **interchangeable cushions and forehead pads**. The cushions are available in various materials and sizes. **The patented adjustable exhalation port** increases comfort by allowing the patient to select the direction the air exits the mask – no more blowing on the patient's chest or bed partner. **The lightweight mask is comfortable and provides a good seal.**



DeVilbiss Serenity  
CPAP Mask  
9352D

### Performance Features

- Patented adjustable forehead rest takes pressure off the bridge of the nose for maximum patient comfort
- Patented adjustable vent allows patient to direct escaping air away from their chest or bed partner
- Snap together components are easy to disassemble and reassemble and minimizes likelihood of breakage
- Solid and reliable thermoplastic frame will hold up to repeated cleanings



## Design Features

- 100% latex free
- Soft replaceable cushion provides a comfortable low pressure seal. Available in standard and shallow sizes and silicone or gel material
- Soft replaceable silicone or gel forehead pads
- Tubing runs over the head for improved tubing management and patient comfort
- Minimal headgear offers a secure fit and features a quick release buckle
- Supplemental oxygen port

## Reimbursement information

Product Description	Code Description	HCPCS Code	Frequency	Allowable Range
Serenity Mask	Nasal Application Device (mask or cannula type), used with positive airway pressure device, with or without head strap	A7034	1 per 3 months	\$99.00 - \$117.64 initial set-up
Replacement Cushion	Replacement cushion for Nasal Application Device, each	A7032	2 per 1 month	\$34.45 - \$40.53
Replacement Headgear	Headgear	A7035	1 per 6 months	\$33.79 - \$39.75
Chin Strap	Chinstrap	A7036	1 per 6 months	\$15.47 - \$18.20

## Accessory and Replacement Parts

DEVILBISS PART #	PRODUCT DESCRIPTION
9352D	Serenity Mask with Standard Size Silicone Cushion
9352G	Serenity Mask with Standard Size Gel Cushion
9352S	Serenity Mask with Shallow Size Silicone Cushion
9352GS	Serenity Mask with Shallow Size Gel Cushion
9352D-601	Replacement Standard Size Silicone Cushion
9352G-601	Replacement Standard Size Gel Cushion
9352S-601	Replacement Shallow Size Silicone Cushion
9352GS-601	Replacement Shallow Size Gel Cushion
9352D-604	Silicone Forehead Pads
9352G-604	Gel Forehead Pads
9352D-619	Headgear



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501 USA  
 800.338.1988 | 814.443.4881 | [www.DeVilbissHealthcare.com](http://www.DeVilbissHealthcare.com)