

t h e c o m f o r t s e r i e s ™

 **Direct Home  
MEDICAL**

Presented With Compliments From  
[www.DirectHomeMedical.com](http://www.DirectHomeMedical.com)



The bridge of the nose.  
Okay, enough of that sore subject.



To find a comfortable nasal mask just follow your nose. That's where you'll see Respironics Simplicity™, a part of the Comfort Series of masks. And "see" is the operative word. It's small. Lightweight. Worn just on the tip of the nose, it eliminates pressure on the bridge of the nose. And patients can wear it and their glasses at the same time. They may find that the most comfortable feature of all.

## Respironics® Simplicity™ Nasal Mask: Big Comfort in a Small Package.

With its innovative design and streamlined construction, Simplicity is ideal for use in the treatment of adult Obstructive Sleep Apnea and for noninvasive ventilation in the home or hospital. It's a perfect fit with Respironics Power Programs™ for Providers, helping to ensure long-term compliance with CPAP and bi-level positive pressure therapy.

Simplicity's small, lightweight cushion does not contact the bridge of the nose, eliminating pressure and irritation. This allows patients to wear eyeglasses and eliminates airflow to the eyes – a contributing factor in contracting conjunctivitis.

- The unique built-in exhalation port produces quiet exhalation flow for enhanced comfort and undisturbed sleep for the patient's bed partner.
- Two sizes fit a wide range of patients for easy mask sizing and inventory management. A unique one-size-fits-most headgear is packaged with the mask.
- The preassembled mask simplifies use and reduces "lost" parts. The mask assembly and headgear are the only pieces that need to be cleaned.

Ordering Information	Part Number
Small Mask with Headgear	1002757
Medium Mask with Headgear	1002759
Replacement Parts	
Headgear	1002764
Small Replacement Cushion	1003658
Medium Replacement Cushion	1003659
Tubing Assembly	1004209

The **exhalation port** results in quiet air flow, and its innovative built-in design eliminates the need for additional valves or ports.

**360° swivels** at the top and bottom of the tubing allow the tubing to twist more freely as the patient moves.

The **mask cushion** is made of soft, pliable silicone that gently contours around the tip of the patient's nose to provide a snug seal.

**Flexible, lightweight tubing** is permanently attached to the mask and connects easily to standard circuit tubing.

The **Stability Clip** on the tubing may be attached to sleepwear or to the headgear for added support. A Velcro® tab secures placement of the Stability Clip on the tubing.

The **Simplicity Nasal Mask Headgear** is specially designed for a secure, comfortable fit. Made of soft, flexible neoprene, the Headgear features a unique five-strap design that evenly distributes tension, adds stability, and helps to prevent mask leaks. Velcro® tabs allow easy adjustment for a custom fit.



Customer Service: 1-800-345-6443 • 724-387-4000  
Respironics Europe: +33-1-55-60-19-80  
Respironics Asia Pacific: +81-3-5280-9611

Simplicity, Comfort Series and Power Programs are trademarks of Respironics, Inc.  
Velcro® is a trademark of Velcro Industries.  
Latex-Free • Patent Pending  
©2004 Respironics, Inc.

**CAUTION:** U. S. federal law restricts this device to sale by or on the order of a physician.

1025656 BD 12/16/04 MCI 4100117

