

Swift LT Nasal Pillows System

Designed for your lifestyle, the Swift LT nasal pillows system is light on the face, easy to use and fit, and provides personalized comfort and unparalleled quietness for you and your bed partner.



RESMED

DirectHome
MEDICAL

Presented With Compliments From
www.DirectHomeMedical.com

MyRESMED

Background | Products | Resources | Therapy Tips | My Therapy

Sign In

Products & Services

MyResMed now features detailed information about all of our products! [Learn more...](#)

- New! The S8/Swift system
- New! The ResMed pediatric system
- Browse product catalog

Therapy Tips

Get the most out of your therapy - take advantage of these useful tips and tricks. [Learn more...](#)

- Choosing the right mask
- Choosing the right PAP
- Traveling with PAP

Mask Replacement

Did you know your insurance may cover a new mask every 6 months? [Learn more...](#)


- Free replacement
- Mask reminders
- Warranty registration

ResMed Corp Poway, CA, USA +1 858 746 2400 or 1 800 424 0737 (toll free), ResMed Ltd Bella Vista, NSW, Australia +61 (2) 8884 1000 or 1 800 658 189 (toll free). Offices in Austria, Belgium, Brazil, China, Finland, France, Germany, Greece, Hong Kong, India, Ireland, Italy, Japan, Malaysia, Mexico, Netherlands, New Zealand, Norway, Portugal, Singapore, South Africa, Spain, Sweden, Switzerland, United Kingdom (see website for details). Swift is a trademark of ResMed Ltd. ©2008 ResMed. 1012407/1 08 04

RESMED

Swift™ LT NASAL PILLOWS SYSTEM Light touch, Easy fit





Featuring...

Light touch, quiet comfort

- **Compact and streamlined design**
Minimizes facial contact while offering stability with a secure seal (ideal for active sleepers)
- **Soft, flexible and stable headgear**
Flexes and adjusts to your facial width and minimizes facial marks
- **Dual-wall nasal pillows**
Sit on a flexible pillow base and seal softly to reduce airflow into the nasal passages for maximum comfort
- **Lightweight design**
Weighs only 2.3 oz (67 g), one of the lightest on the market
- **Unique vent design**
Delivers super quiet performance, making Swift™ LT the quietest on the market



Personalized fit and comfort

- **Rotating barrel**
Allows you to personalize the most comfortable fit the first time, every time for your nasal pillows without over-tightening your headgear
- **Flexible lightweight tubing**
360° swivel enables multiple sleeping positions, including side or back
- **Optional tube retainer**
Clicks onto top of headgear for multiple tubing positions including side and back sleeping



Easy to fit, easy to use

- **First-time fit**
Ready to fit straight from package with minimal adjustments
- **Three pillow sizes (S, M, L) included**
Allows you to find the perfect fit
- **Fewer parts**
Allows for trouble-free cleaning and assembly
- **Convenient drawstring bag**
Fits your lifestyle and travel needs

