


INOGEN ONE[®] G4 System



Inogen One[®]
AC Power Supply*
Item # BA-401



Inogen One[®] G4
DC Power Cable*
Item # BA-306



Inogen One[®] G4
Lithium Ion Battery*
Up to 2.7 hours run time
Item #BA-400



Inogen One[®] G4
Carry Strap*
Item # CA-401



Inogen One[®] G4
Carry Bag*
Item # CA-400



Inogen One[®] G4
Extended Life Lithium Ion Battery**
Up to 5 hours run time
Item #BA-408

The Inogen One[®] G4 System (IS-400):

- Compact and Lightweight
- Noticeably Quiet Operation
- Intelligent Oxygen Delivery
- Designed for 24/7 Use

*Included with system

**Accessory sold separately

**Only
2.8 pounds!**



Oxygen. Anytime. Anywhere[®]

The Inogen One G4 delivers the independence of a portable oxygen concentrator in one of the smallest, lightest, and quietest packages available to the oxygen user today. With flow settings from 1-3, the Inogen One G4 is designed to increase independence for most users 24/7. It can be charged at home or on the go - giving you the freedom of Oxygen. Anytime. Anywhere. [®]

The Inogen One[®] G4 is your one solution, your single solution, for oxygen at home or away.

| | |
|--------------------|--|
| Weight | 2.8 lbs. (includes single battery) |
| Size | Length: 5.91" Width: 2.68" Height: 7.2" (includes single battery) |
| Oxygen Flow | Pulse dose delivery system Three flow settings 1, 2, 3 |
| Power | AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car |
| Battery | Duration (single battery): Up to 2.7 hours Recharge approximately 3 hours with AC or DC power Duration (double battery): Up to 5 hours Recharge approximately 5 hours with AC or DC power |
| Noise Level | 40 dBA* |
| Warranty | 3 year standard warranty** |
| Operation | Simple control functions and easy-to-read LCD display |
| Use | Designed for 24/7 use and conforms to all applicable FAA requirements for POC use on board an aircraft |

*At flow setting 2

**See Inogen One[®] Limited Warranty Statement for complete warranty information

Manufactured by:



326 Bollay Drive
Goleta, California 93117
www.inogen.com
www.inogen.com/hme-providers

