

# Sleeping Is a Breeze With Puritan Bennett® CPAP Systems.



Puritan Bennett has been a leading provider of CPAP systems for over 15 years. The *GoodKnight*® family of CPAP systems is known for its reliability and whisper quiet operation – letting sleep apnea patients get the rest they need.

The *GoodKnight* 420 Series is the world's smallest and lightest CPAP device. *GoodKnight* 420 units have special travel-friendly features that make it easy for patients to comply with therapy wherever they go:

- > Automatic voltage regulator
- > Optional D/C power supply
- > Altitude compensation feature

GoodKnight Series Models Include:

GoodKnight 420G  
GoodKnight 420S  
GoodKnight 420E

Features	420G Entry Level	420S W/Compliance Tracking	420E Auto-CPAP
Compliance and Hour Meters	+	+	+
Pressure Sensor	+	+	+
Automatic Leak Compensation	+	+	+
Automatic Altitude Compensation	+	+	+
Ramp Feature	+	+	+
Compliance Data Can Be Downloaded		+	+
SilverLining® Data Management Software Option		+	+
Detects Apnea, Hypopnea and Snoring Events		+	+
Reacts to Apnea, Hypopnea and Snoring Events			+
Detects and Reacts to Inspiratory Flow Limitation			+
Detects Apnea with and without Cardiac Oscillations			+

## Breeze SleepGear™

CPAP Interface System

Ordering Information	
Description	Order Number
<b>Breeze SleepGear Interface Systems With Nasal Pillows*</b>	
<i>Breeze SleepGear</i> Interface System with Medium and Large Nasal Pillows	Y-101400-00
<i>Breeze SleepGear</i> Interface System with Small Nasal Pillows	Y-101400-S
<i>Breeze SleepGear</i> Interface System with Medium Nasal Pillows	Y-101400-M
<i>Breeze SleepGear</i> Interface System with Large Nasal Pillows	Y-101400-L
<i>Breeze SleepGear</i> Interface System with Extra Large Nasal Pillows	Y-101400-XL
<i>Breeze SleepGear</i> Interface System with Small Dilators	Y-101400-SD
<i>Breeze SleepGear</i> Interface System with Medium Dilators	Y-101400-MD
<i>Breeze SleepGear</i> Interface System with Large Dilators	Y-101400-LD
<b>Breeze SleepGear Interface Systems With DreamSeal Mask*</b>	
<i>Breeze SleepGear</i> Interface System with Standard <i>DreamSeal</i> Mask	Y-102616-00
<i>Breeze SleepGear</i> Interface System with Large <i>DreamSeal</i> Mask	Y-103059-00A
<i>Breeze SleepGear</i> Interface System with Shallow <i>DreamSeal</i> Mask	Y-103074-00A
* Interface Systems include halo straps, 15-mm tubing and airway assembly	
<b>Nasal Pillows Accessories/Replacement Parts</b>	
Nasal Pillows Shell with pressure tap (for sleep lab use)	Y-101987-00
Nasal Pillows, Small, blue, 1 pair	S-616323-00A
Nasal Pillows, Medium, clear, 1 pair	S-616324-00C
Nasal Pillows, Large, green, 1 pair	S-616325-00A
Nasal Pillows, Extra Large, lavender, 1 pair	Y-101969-00
Dilator, Small, gray, 1 pair	S-616355-00C
Dilator, Medium, white, 1 pair	S-616356-00C
Dilator, Large, yellow, 1 pair	S-616357-00C
Nasal Airway Assembly without pillows, includes shell w/vent, 15-mm tubing, swivel adapter and hose guide	Y-102444-00
Nasal Pillows Shell, includes vent	Y-101699-00
Nasal Pillows Exhalation Vent, 3/pkg	Y-101974-00
Pillow Sizing Gauge, 10/pkg	S-616468-PG
<b>DreamSeal Mask Accessories/Replacement Parts</b>	
<i>DreamSeal</i> Mask Assembly, includes shell with vent, 15-mm tubing, swivel adapter, hose guide	Y-102620-00
<i>DreamSeal</i> Mask Assembly with seal, Standard	Y-102617-00
<i>DreamSeal</i> Mask Assembly with seal, Large	Y-103061-00A
<i>DreamSeal</i> Mask Assembly with seal, Shallow	Y-103076-00A
<i>DreamSeal</i> Shell, includes vent	Y-102622-00
<i>DreamSeal</i> Shell with pressure tap (for sleep lab use)	Y-102637-00
<i>DreamSeal</i> Exhalation Vent, 3/pkg	Y-102638-00
<i>DreamSeal</i> Mask, Standard	Y-102639-00
<i>DreamSeal</i> Mask, Large	Y-103060-00A
<i>DreamSeal</i> Mask, Shallow	Y-103075-00A
<b>Breeze SleepGear Universal Components</b>	
<i>Breeze SleepGear</i> Assembly, includes spine, cradle and outrigger	Y-102640-00
<i>Breeze SleepGear</i> Cradle Assembly	Y-102950-00
Halo Straps, 1 pair	Y-101975-00
15-mm Tubing and 15-mm/22-mm Swivel Adapter	Y-101976-00
Hose Guide Kit, includes hose guide, clip and halo straps	Y-102833-00
Swivel Adapter	S-616329-00B

**tyco**  
Healthcare

**Puritan Bennett**

4280 Hacienda Drive  
Pleasanton, CA 94588  
Tel 925.463.4000  
Toll Free 1.800.635.5267  
www.puritanbennett.com

Tyco Healthcare UK LTD.  
154 Fareham Road  
Gosport, UK PO13 0AS  
Tel +44.1329.224000

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician. Specifications are subject to change without notice. Patent number 4782835.  
© 2004 Nellcor Puritan Bennett Inc. All rights reserved. A.d 0530v6-0504

Specifications
<b>Materials of Construction</b>
Nasal Pillows: Medical grade silicone <i>DreamSeal</i> Mask: Perimeter seal is medical grade silicone Straps: Nylon and spandex All <i>Breeze SleepGear</i> components are latex free.
<b>Pressure Range</b>
3-35 cm H <sub>2</sub> O
<b>Hose Attachment</b>
Fits 22 mm OD connection

*Breeze SleepGear* With Nasal Pillows Assembly



*Breeze SleepGear* With *DreamSeal* Assembly



Nasal Pillows

*DreamSeal* Mask



A.b 1751-0504 ST00900

This unique interface can be used with the Nasal Pillows or *DreamSeal*® Mask assembly.

PURITAN BENNETT

# Breeze SleepGear™

CPAP Interface System



Presented With Compliments From  
www.DirectHomeMedical.com

Breathe easy. Sleep better.

# Breeze SleepGear™ Interface Your First Choice in Fit and Comfort.

## Breeze SleepGear CPAP Interface System

The unique *Breeze SleepGear*™ CPAP Interface System from Puritan Bennett sets new standards in comfort and versatility for sleep apnea patients. Unlike traditional headgear and masks, there are no complicated or confining straps to contend with. The streamlined design contours the head and floats over the face, giving patients the freedom to sleep in virtually any position.

### Two-in-One Versatility

The innovative *Breeze SleepGear* system lets patients easily switch between the Puritan Bennett Nasal Pillows and *DreamSeal*® Mask, so they aren't limited to the same type of interface night after night. The patented Nasal Pillows are available in seven sizes, offering a wide choice in fit. The *DreamSeal* Mask comes in three sizes.

PURITAN BENNETT

For more information call 1.800.635.5267  
or visit [www.puritanbennett.com](http://www.puritanbennett.com)



## Fitting the Breeze SleepGear Interface with Nasal Pillows



## For Hard-to-Fit Patients and Every Patient in Between.

Make the *Breeze SleepGear* system your first choice in headgear for its ability to accommodate a wide range of patients. Simply slide the Nasal Pillows or *DreamSeal* Mask assembly into the *Breeze SleepGear* interface, and take advantage of the key adjustment points shown below.

- > New clip (Nasal Pillows only) for patients who want more stability.
- > Multiple adjustments and size options make it easy to get just the right fit.
- > Versatile headgear works with the Nasal Pillows or *DreamSeal* Mask.
- > No confining, hard-to-adjust straps.
- > Ultra-quiet exhalation vent.
- > Patients can sleep in virtually any position.



# New & Improved Features:

For over five years, users of the *Breeze SleepGear* CPAP Interface System have shared their gratitude for its comfort and performance. They also shared feedback and suggestions that we listened to in adding new features to make *Breeze SleepGear* system even better.

- > **Durability**—The hose guide for the Nasal Pillows assembly has an improved design and features a stronger material to ensure greater longevity.
- > **Fit**—A new clip has been added to the Nasal Pillows assembly to give patients greater flexibility in positioning the halo straps for a secure fit.
- > **Comfort**—The cradle pad has thicker padding to increase comfort for a good night's sleep.
- > **Convenience**—A cuff has been added to the swivel adapter for improved convenience.

# Sleep Lab Certified

The *Breeze SleepGear* system offers convenience and efficiency in the sleep lab. When the headgear is applied before the patient falls asleep, it is easy to connect the *DreamSeal* Mask or Nasal Pillows assembly for titration with minimal disruption to the patient.

"Having a 'strapless system' for titrations in the labs has reduced the number of lead repairs done during the night, thanks to *Breeze SleepGear*. And CPAP instruction and setups can be done twice as fast due to its simplicity."

— Margaret Smith, RPSGT, CRTT