

Interfaces

FlexiFit™ 431 Full Face Mask



The *FlexiFit™ 431* Full Face Mask with Under Chin Design builds on the foundation of the unique features common to the *FlexiFit™* Series of nasal masks.

The *FlexiFit™* Technology has been incorporated into the contoured Silicone Seal to allow a large range of nasal bridge movement and simple one step auto adjustment fitting. The unique *Glider™* Strap provides freedom of movement while maintaining a seal. The crown strap and Under Chin Design provide extra stability and maintain a fit even with jaw movement. The Easy-Release Strap allows simple removal of the mask and the “Fit and Forget” Headgear ensures no on-going adjustments are required.

The *FlexiFit™ 431* Full Face is supplied fully assembled, with a second size of Silicone Seal in the box.

Product Code	Description	Quantity
HC431	<i>FlexiFit™ 431</i> Full Face Mask includes: <ul style="list-style-type: none">• Silicone Seal (Small)• Silicone Seal (Medium)• Silicone Seal (Large)• Headgear with Easy Release Cord• <i>Glider™</i> Strap	1 1 1 1 1

BREATH-O-PRENE™ is a registered trade mark of AccuMED Technologies, Inc



Presented With Compliments From
www.DirectHomeMedical.com

Fisher & Paykel
HEALTHCARE

FlexiFit™ 431 Mask Spare Parts



Elbow Kit


Product Code	Description	Quantity
400HC203	FlexiFit™ 431 Elbow & Vented Non Re-Breathing Valve	Each
400HC201	FlexiFit™ 431 Elbow & Non Re-Breathing Valve (Aus only)	Each

Warning: Do not use with Non-Invasive Ventilation models.



Silicone Seal

Product Code	Description	Quantity
400HC112	Silicone Seal FlexiFit™ 431 Small	Each
400HC104	Silicone Seal FlexiFit™ 431 Medium	Each
400HC103	Silicone Seal FlexiFit™ 431 Large	Each



Forehead Pads

Product Code	Description	Quantity
900HC423	Forehead Pads - fits FlexiFit™ 405, 407, 431	2 Pack



Oxygen/Pressure Port Connector

Product Code	Description	Quantity
900HC452	Oxygen/Pressure Port Connector	10 Pack



Mask Headgear

Product Code	Description	Quantity
400HC302	FlexiFit™ 431 Headgear with Crown Strap	Each



Mask Kit

Product Code	Description	Quantity
400HC503	FlexiFit™ 431 Mask – No Headgear	Each
400HC518	FlexiFit™ 431 Mask – No Headgear Australia Only	Each