



Presented With Compliments From  
[www.DirectHomeMedical.com](http://www.DirectHomeMedical.com)

# Mirage Quattro™

FULL FACE MASK


## FOURTH GENERATION TECHNOLOGY

**Wide fit range:** fourth generation full face mask from ResMed fits more than 95% of users.\*

**Mirage™ dual-wall cushion with spring air technology:** unique cushion accommodates movement during sleep and provides an effective solution for jaw drop.

**MicroFit dial:** new forehead support adjustment offers personalized fitting and comfort.

**Light, quiet and easy to use:** these benefits are optimized with streamlined mask design, improved vent design and a range of user-friendly features.





## TECHNOLOGY HIGHLIGHT

### MicroFit dial

A choice of 24 positions at your fingertips—providing a comfortable, effective fit for over 95% of users\*

- **Adjustable:** can be adjusted while on the user's face
- **Precise:** each turn of the dial adjusts the forehead support in small increments
- **Progressive:** moves in towards the face, allowing the user to sense when the fit is most comfortable for them
- **Works in combination with the cushion:** as the dial adjusts the angle of the mask to the face, the cushion adapts to enhance the seal where it's needed most—without the need to overtighten the headgear

## USER BENEFITS

### Personalized comfort

- MicroFit dial adjusts the mask to suit user's profile
- Dual-wall cushion with spring air technology adapts to facial contours and accommodates movement during sleep (including jaw drop)

### Light and quiet

- Streamlined design feels light on the face
- Improved vent design reduces noise

### Making it easier

- MicroFit dial and spring air cushion technology help achieve good mask fit—first time, every time
- Fewer parts for easy handling, cleaning, assembly and disassembly

## PROFESSIONAL BENEFITS

### Fit range increased to over 95% of users\*

- MicroFit dial and spring air cushion technology adjust the mask to suit more facial structures
- New XS mask size

### First time fit

- First-time fit reduces callbacks and mask refitting
- Improved cushion technology seals quickly
- Simplified parts for quick and easy assembly and disassembly

### Reduced inventory

- 33% fewer product codes than Ultra Mirage™ Full Face Mask
- Fewer mask parts

### Comprehensive user support package

- Travel and storage bag
- Mask care DVD
- Quick-fitting guide
- User guide



## PRODUCT CODES

### Mirage Quattro Full Face Mask

X-Small	<b>61200</b>
Small	<b>61201</b>
Medium	<b>61202</b>
Large	<b>61203</b>

## REIMBURSEMENT CODES

### Mirage Quattro Full Face Mask

HCPSC Code	Descriptor
<b>A7030</b>	Full face mask used with positive airway pressure (PAP) device, 1 per 3 months
<b>A7031</b>	Replacement cushion for full face mask, 1 per month
<b>A7035</b>	Headgear used with positive airway pressure (PAP) device, 1 per 6 months

\* Source: internal ResMed anthropometric data